

**EUROPEAN AND MEDITERRANEAN PLANT PROTECTION ORGANIZATION**  
**ORGANISATION EUROPEENNE ET MEDITERRANEENNE POUR LA PROTECTION DES PLANTES**  
**Summary sheet of validation data for a diagnostic test**

The EPPO Standard PM 7/98 *Specific requirements for laboratories preparing accreditation for a plant pest diagnostic activity* describes how validation should be conducted. It also includes definitions of performance criteria.

<b>Laboratory contact details</b>	Finnish Food Authority / Plant Pest Section Mustialankatu 3, 00790 Helsinki, Finland
<b>Short description of the test</b>	Detection of 'Candidatus Liberibacter solanacearum' using real-time PCR test
<b>Date, reference of the validation report</b>	2018-04-03 - Validatoin_report_Evira7629, ID DOC-1280-6573-fi
<b>Validation process according to EPPO Standard PM7/98?</b>	yes
<b>Is the lab accredited for this test?</b>	no
<b>Was the validated data generated in the framework of a project?</b>	
<b>Description of the test</b>	
<b>Organism(s)</b>	'Candidatus Liberibacter solanacearum' (LIBEPS)
<b>Detection / identification</b>	detection
<b>Matrix(ces) tested</b>	Tubers Tuber
<b>Plant species tested</b>	Solanum tuberosum
<b>Method(s)</b>	Extraction Molecular Extraction DNA RNA Molecular real time PCR
<b>Method: Extraction</b>	
<b>Reference of the test description</b>	
<b>Other information</b>	
<b>Other details on the test</b>	Potato tuber sample can be taken from a single tuber, or from the heel end core sample for ring rots (LAB4001) directly or after the ring rots have been extracted. Heel end cores are homogenized in Bioreba extraction bag without the buffer, 25 ml PBST is added and the bacteria are extracted by rubbing the buffer into the tuber pulp inside the bag. The bacteria are concentrated from the extract by 2-step centrifugation. In the first step (clarification) the coarse plant material is removed from the sample by 100×g 5 min centrifugation at 1 / 4 10 °C, the bacteria are then pelleted from the supernatant (~8.5 ml) by 7000×g 5 min centrifugation at 10 °C. Bacterial pellet is suspended in 500 µl PBST buffer, 200 µl of the

suspension is used for DNA extraction and the rest can be frozen in case the analysis needs to be repeated. A single tuber sample is taken from the heel end or bud site, homogenized with Homex 6 and extracted 1:2 to PBST buffer. It is important to get as much vascular tissue as possible into the sample. Symptomatic sample does not need concentration. Carrot petioles and potato stolons are homogenized in in Bioreba extraction bag with Homex 6 and extracted 1:2 in PBST. 1 cm pieces from the petioles are cut circa 3 cm above the root. 200 µl of the unconcentrated extract is used for DNA extraction. Carrot seeds are washed before grinding to remove the seed dressing, in addition the washing softens the seeds and thus aids the grinding. One gram of seeds is weighed into 50 ml centrifuge tube and 30–40 ml 0.5 % Triton X-100 is added. Larger volume reduces the flow of liquid in the mixing and results in poor washing result. The tube is shaken 30 minutes on a platform rocker. The liquid is carefully poured out, the tube is filled with UP-water, shaken vigorously and the water is carefully decanted. The rinsing is performed three times total, on the third time the seeds are poured with the water into Bioreba extraction bag and the water is subsequently decanted from the bag. The seeds are ground to fine powder with Homex 6, 10 ml PBST is added and the bacteria are extracted by rubbing the buffer into the seed powder inside the bag. 200 µl of the extract is used for DNA extraction. Color from seed dressing may remain in the sample, it does not interfere with the extraction.

**Method: Molecular Extraction DNA RNA**

**Reference of the test description**

**Kit**

**Is a kit used**

yes

**Manufacturer name**

BIONOBILE

**Specify the kit used**

QuickPick Plant DNA kit

Kit used following the manufacturer's instructions?

**Other information**

**Other details on the test**

DNA extraction with KingFisher Flex and QuickPick™ Plant DNA-kit (BioNobile) 1. Add 300 µl Lysis Buffer and 20 µl proteinase K to 200 µl sample 2. Incubate 30 minutes in thermal shaker 700 rpm at 65°C 3. Pipet KingFisher Deepwell plates (Cat. 95040450) during incubation Sample: 500 µl Binding Buffer + 20 µl Magnetic Particles Wash 1: 900 µl Wash Buffer Wash 2\_1: 800 µl Wash Buffer Wash 2\_2: 800 µl Wash Buffer Elution: 50 µl Elution Buffer 4. Centrifuge the samples 5 minutes 18000xg at room temperature. Pipet 450 µl of the supernatant carefully to the Sample plate wells avoiding the sediment. 5. Start KingFisher Flex and

	run PURE_DNAPlant_Flex96 program 2 / 4 6. Pipet the DNA from Elution plate to eppendorf tube and use for real-time PCR analysis. The DNA can be kept short time (1-2 weeks) in the fridge or longer times in the freezer for further use.
<b>Method: Molecular real time PCR</b>	
<b>Reference of the test description</b>	
<b>As or adapted from an EPPO diagnostic protocol</b>	yes
<b>EPPO Diagnostic Protocol name</b>	PM 7/143 ' <i>Candidatus Liberibacter solanacearum</i> ' (version 1)
<b>As or adapted from an IPPC diagnostic protocol</b>	no
<b>Is the test modified compared to the reference test</b>	yes The primers and probes adopted from Li et al. (2009) Journal of Microbiological Methods 78:59-65 (universal HLBr primer and HLbP probe, CLso specific LsoF primer) or Teresani et al. (2014) Phytopathology 104: 804-811 (universal CaLsppF and CaLsppR primers, CLso specific CaLsolP probe) were used in conjunction with in-house developed DNA extraction method and modified real-time PCR (qPCR). Deviations from the references: - Sample preparation - DNA extraction - PCR reagents and reaction volume - Primer and probe concentrations - BHQ1-quencher is used instead of TAMRA for CaLsolP probe - PCR program (Teresani et al. 2014 is also used for Li et al. 2009 primers)
<b>Other information</b>	
<b>Other details on the test</b>	PCR reaction for one well Reagent Conc. Volume (µl) Final conc. Nuclease free water 4.0 Maxima Probe qPCR Master Mix 2× 12.5 1× CaLsppF OR LsoF OR COX-F 5 µM 1.5 0.3 µM CaLsppR OR HLBr OR COX-RW 5 µM 1.5 0.3 µM CaLsolP OR HLbP OR COX-P 5 µM 0.5 0.1 µM Total 20.0 DNA 5.0 Total 25.0 PCR An initial step at 95°C for 10 min followed by 45 cycles (95°C for 15 s, and 60°C for 1 min; Teresani et al. 2014).
<b>Are the performance characteristics included in the EPPO diagnostic protocol?</b>	<b>no</b>
<b>Performance Criteria :</b>	
<b>Organism 1.:</b>	<b>'Candidatus Liberibacter solanacearum'(LIBEPS)</b>
<b>Analytical sensitivity</b>	
<b>What is the smallest amount of target that can be detected reliably?</b>	0.01% contamination rate or 10 target copies per PCR reaction
<b>Diagnostic sensitivity</b>	
<b>Proportion of infected/infested samples tested positive compared to results from the standard test, see appendix 2 of PM 7/98</b>	N.A.

<b>Standard test(s)</b>	N.A.
<b>Analytical specificity - inclusivity</b>	
<b>Number of strains/populations of target organisms tested</b>	4: CLso-positive <i>Daucus carota</i> subsp. <i>sativus</i> midrib (83-2013, haplotype C, Minna Haapalainen, University of Helsinki) and seeds [cv. Maestro (11-2016, haplotype D, Mikko Lehtonen, Evira)], <i>Solanum tuberosum</i> tuber frozen (2-2011, haplotype B, Anne Nissinen, LUKE) and fresh material [cv. Atlantic (1-7-2017, haplotype B, Rodney Cooper, USDA)]
<b>Specificity value</b>	100%
<b>Analytical specificity - exclusivity</b>	
<b>Number of non-target organisms tested</b>	4: <i>Ralstonia solanacearum</i> NCPPB 4156, <i>Clavibacter michiganensis</i> var. <i>sepedonicus</i> NCPPB 4053, and <i>Dickeya solani</i> S0432-1 and <i>Pectobacterium carotovorum</i> SCC1 (Minna Pirhonen, UH).
<b>Specificity value</b>	100%
<b>Diagnostic Specificity</b>	
<b>Proportion of uninfected/uninfested samples (true negatives) testing negative compared to results from a standard test</b>	N.A.
<b>Specify the test(s)</b>	N.A.
<b>Reproducibility</b>	
<b>Provide the calculated % of agreement for a given level of the pest (see PM 7/98)</b>	100% for HLBp, 66% for CaLsoLP
<b>Repeatability</b>	
<b>Provide the calculated % of agreement for a given level of the pest (see PM 7/98)</b>	100%
<b>Test performance study</b>	
<b>Test performance study?</b>	no
The following complementary files are available online:	
	<ul style="list-style-type: none"> <li>• <a href="#">SOP english</a></li> <li>• <a href="#">Validation report</a></li> </ul>

Creation date: 2018-11-07 00:00:00 - Last update: 2021-05-14 17:23:53